

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB09186 - Goat Anti-HIP1R/ HIP12 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: HIP1R, huntingtin interacting protein 1 related, FLJ14000, HIP12, HIP3, ILWEQ, KIAA0655, MGC47513, huntingtin interacting protein 12, huntingtin

interacting protein-1-related **Official Symbol:** HIP1R

Accession Number(s): NP_003950.1

Human GeneID(s): 9026

Immunogen

Peptide with sequence C-RQDHQLDKKDGIYP, from the C Terminus of the protein sequence according to NP_003950.1.

Please note the <u>peptide</u> is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Dog, Cow, Pig