



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05524 - Goat Anti-HMG20B (BRAF35) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: HMG20B, high-mobility group 20B, SOXL, HMGX2, BRAF25, BRAF35, SMARCE1r, Sox-like transcriptional factor, SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily E, member 1-related

Official Symbol: HMG20B

Accession Number(s): NP_006330

Human GeneID(s): [10362](#)

Immunogen

Peptide with sequence C-QDPPNYSLRGAVP, from the C Terminus of the protein sequence according to NP_006330.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 35-40kDa band seen in Human Heart lysate (Predicted MW of 38kDa according to NP_006330). Recommended for use at 2-4µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human