



## EB2ND-50-3001-HRP - Goat anti-Human IgG antibody HRP Conjugated

Size: 500µg



### UK Office

#### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

### Immunogen

Purified whole human IgG (Horseradish Peroxidase/ HRP conjugated).

### Purification

Ammonium sulphate precipitation of serum followed by antigen affinity purification.

### Buffer

Phosphate buffer saline, pH 7.2, with 10 mg/ml BSA and 0.01% Gentamicin sulphate as preservative. Quantity: 0.5ml.

### Storage and Handling

Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles.

### Applications

Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based immunoassays. Recommended dilutions:

ELISA: 1:50000 WB: 1:25000-1:50000. Please optimize dilutions for your particular assay conditions. Starting test dilution: 0.5µg/ml for ELISA and WB.