



EB2ND-50-3001 - Goat anti-Human IgG antibody

Size: 500µg



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre

77 Heyford Park

Upper Heyford

Oxfordshire

OX25 5HD

UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

Immunogen

Purified whole human IgG (unconjugated).

Purification

Ammonium sulphate precipitation of serum followed by antigen affinity purification.

Buffer

Phosphate buffer saline, pH 7.2, with 0.02% Sodium azide as preservative. Quantity:
0.5mg in 0.5mg/ml.

Storage and Handling

Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize
freeze thaw cycles.

Applications

Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based
immunoassays. Recommended dilutions:

ELISA: 1:128000 WB: 1:3000 to 1:10000. Please optimize dilutions for your particular
assay conditions. Starting test dilution: 0.5µg/ml for ELISA and WB.