



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB2ND-50-3001 - Goat anti-Human IgG antibody

Size: 500µg



### Immunogen

Purified whole human IgG (unconjugated).

### Purification

Ammonium sulphate precipitation of serum followed by antigen affinity purification.

### Buffer

Phosphate buffer saline, pH 7.2, with 0.02% Sodium azide as preservative. Quantity:  
0.5mg in 0.5mg/ml.

### Storage and Handling

Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize  
freeze thaw cycles.

### Applications

Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based  
immunoassays. Recommended dilutions:

ELISA: 1:128000 WB: 1:3000 to 1:10000. Please optimize dilutions for your particular  
assay conditions. Starting test dilution: 0.5µg/ml for ELISA and WB.