



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB2ND-50-3002 - Goat anti-human IgM antibody

Size: 500ug

Immunogen

a 42-aa long peptide from the heavy chain of human Immunoglobulin M (unconjugated).

Purification

Ammonium sulphate precipitation of serum followed by antigen affinity purification.

Buffer

10 mMTris, PH 8.2, with 150 mMNaCl, 10 mg/ml BSA and 0.1% sodium azide

Storage and Handling

Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles.

Applications

Protein ELISA: Human serum IgM detection at 1: 200 dilution.

Peptide ELISA: Antibody detection at 1:128000 dilution.

Immunofluorescence: Immunofluorescence tested on Daudi cells. Recommended dilution: 1:50 to 1:100.

Western blotting : Approx 75 kDaband observed human plasma and serum.

Recommended concentration: 0.1-1µg/ml.

Please optimize dilutions for your particular assay conditions. Starting test dilution:1:500 for indirect immunofluorescence.