

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06361 - Goat Anti-HUS1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: HUS1, HUS1 checkpoint homolog (S. pombe), HUS1 (S. pombe) checkpoint homolog

Official Symbol: HUS1

Accession Number(s): NP_004498

Human GeneID(s): [3364](#)

Immunogen

Peptide with sequence KFRKIVDGACLN, from the N Terminus of the protein sequence according to NP_004498.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Western Blot: Preliminary experiments gave no signal but low background in Human Liver and Heart lysates at up to 1µg/ml (note whole cell RIPA lysates not nuclear lysates). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human