



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05282 - Goat Anti-ID2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: ID2, inhibitor of DNA binding 2, dominant negative helix-loop-helix protein, Inhibitor of DNA binding 2, dominant negative

Official Symbol: ID2

Accession Number(s): NM_002166

Human GeneID(s): [3398](#)

Immunogen

Peptide with sequence SELMSNDSKALCG, from the C Terminus of the protein sequence according to NM_002166.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 30kDa band seen in Human Brain lysate.

Recommended for use at 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human