

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB07606 - Goat Anti-IER5 Antibody

Size: 100µg specific antibody in 200µl



**Target Protein** 

Principal Names: IER5, immediate early response 5, MGC102760, SBBI48

Official Symbol: IER5

Accession Number(s): NP\_057629.1

Human GeneID(s): 51278

## **Immunogen**

Peptide with sequence C-QPPSGGEDDDAEE, from the internal region of the protein sequence according to NP\_057629.1.

Please note the peptide is available for sale.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

#### **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:2000.

Western blot: Customer's experiments gave a band at approx 45kDa in lysates of transfected H1299 overexpressing IER5 but not when transfected with empty vector (calculated MW of 33.7kDa according to NP\_057629.1). Recommended concentration: 1-3μg/ml.

## **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human