

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB08349 - Goat Anti-IFNAR2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: IFNAR2, interferon (alpha, beta and omega) receptor 2, IFN-R,

IFN-alpha-REC, IFNABR, IFNARB

Official Symbol: IFNAR2

Accession Number(s): NP_000865.2; NP_997467.1; NP_997468.1

Human GeneID(s): 3455

Important Comments: This antibody is expected to recognise all three reported isoforms

(NP_000865.2; NP_997467.1; NP_997468.1)

Immunogen

Peptide with sequence SYDSPDYTDESCT, from the internal region of the protein sequence according to NP_000865.2; NP_997467.1; NP_997468.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave an approx 28kDa band in Human Uterus lysates after 0.2μg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 37.4kDa according to NP_000865.2. The 28kDa band was successfully blocked by incubation with the immunizing peptide.


IHC: Paraffin embedded Human Prostate and Placenta. Recommended concentration: 5ug/ml.

Flow Cytometry: Flow cytometric analysis of K562 cells. Recommended concentration: 10ug/ml.


Species Reactivity

Tested: Human


Expected from sequence similarity: Human



EB08349 ($5\mu g/ml$) staining of paraffin embedded Human Prostate. Steamed antigen retrieval with citrate buffer Ph 6, AP-staining.



EB08349 ($5\mu g/ml$) staining of paraffin embedded Human Placenta. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.



EB08349 Flow cytometric analysis of paraformaldehyde fixed K562 cells (blue line), permeabilized with 0.5% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (1ug/ml). IgG control:

Unimmunized goat IgG (black line) followed by Alexa Fluor 488 secondary antibody.