



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09691 - Goat Anti-INO80 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: INO80, INO80 homolog (S. cerevisiae), INO80A, INOC1, hINO80, INO80 complex homolog 1, INO80 complex subunit A, homolog of yeast INO80

Official Symbol: INO80

Accession Number(s): NP_060023.1

Human GeneID(s): [54617](#)

Non-Human GeneID(s): 296084 (rat)

Immunogen

Peptide with sequence C-TQEDYDSNHFKAQ, from the internal region of the protein sequence according to NP_060023.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Western blot: Preliminary experiments in Human Brain (Cerebellum) and Human Testis lysates and in lysates of cell lines HeLa and K562 gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Rat Dog, Cow