

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

## EB09140 - Goat Anti-JAM2 / JAM-B Antibody

Size: 100µg specific antibody in 200µl



### **Target Protein**

**Principal Names:** JAM2, junctional adhesion molecule 2, C21orf43, CD322, JAM-B, JAMB, PRO245, VE-JAM, VEJAM, JAM-IT/VE-JAM, OTTHUMP00000096100, junctional

adhesion molecule B, vascular endothelial junction-associated molecule

Official Symbol: JAM2

Accession Number(s): NP\_067042.1

Human GenelD(s): 58494

### **Immunogen**

Peptide with sequence C-NPRLGSQSTNSSYT , from the internal region of the protein sequence according to NP\_067042.1.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments gave an approx 27kDa band in Human Heart, Liver and Placenta lysates after 0.2μg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 33.2kDa according to NP\_067042.1. The 27kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

### **Species Reactivity**

Tested:

Expected from sequence similarity: Human