



UK Office

Everest Biotech Ltd
Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05161 - Goat Anti-JIK: STE20-like kinase Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: JIK, STE20-like kinase, DPK, STE2-like kinase

Official Symbol: TAOK3

Accession Number(s): NP_057365

Human GeneID(s): [51347](#)

Immunogen

Peptide with sequence C-LVTLDFPKEDYR, from the C Terminus of the protein sequence according to NP_057365.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 95 kDa band seen in NIH-3T3 Lysate.

Recommended for use at 1-3 µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human