



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

## US Office

### Everest Biotech c/o Abcore

405 Maple Street, Suite A106  
Ramona,  
CA 92065  
USA

Inquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[usasales@everestbiotech.com](mailto:usasales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB07662 - Goat Anti-KAP1 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** TRIM28, tripartite motif-containing 28, FLJ29029, KAP1, RNF96, TF1B, TIF1B, KRAB-associated protein 1, nuclear corepressor KAP-1, transcriptional intermediary factor 1-beta, tripartite motif-containing 28 protein

**Official Symbol:** TRIM28

**Accession Number(s):** NP\_005753.1

**Human GeneID(s):** [10155](#)

### Immunogen

Peptide with sequence C-KQFNKLTEDKADVQ, from the internal region of the protein sequence according to NP\_005753.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:2000.

**Western blot:** Approx 100kDa band observed in nuclear lysates of cell lines HeLa, Jurkat and K562 (calculated MW of 88.5kDa according to NP\_005753.1). This molecular weight is routinely observed by other sources. Recommended concentration: 1-3µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Mouse, Rat, Dog



EB07662 (1 $\mu$ g/ml) staining HeLa nuclear lysate (35 $\mu$ g protein in RIPA buffer). Primary incubation was 1 hour.  
Detected by chemiluminescence.