



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05470 - Goat Anti-KLF7 (UKLF) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: KLF7, Kruppel-like factor 7 (ubiquitous), UKLF, ubiquitous Kruppel-like transcription factor

Official Symbol: KLF7

Accession Number(s): NP_003700

Human GeneID(s): [8609](#)

Immunogen

Peptide with sequence DVLASYISIFQEL-C, from the N Terminus of the protein sequence according to NP_003700.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 35kDa band seen in Mouse Kidney lysate (predicted MW 33kDa). Recommended for use at 2-4µg/ml.

Species Reactivity

Tested: Mouse

Expected from sequence similarity: Human, Mouse