



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06233-B - Goat Anti-KPNA2 / IPOA1, Biotinylated Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: KPNA2, karyopherin alpha 2 (RAG cohort 1, importin alpha 1), IPOA1, QIP2, RCH1, SRP1alpha, RAG cohort 1, RAG cohort protein 1, SRP1-alpha, importin alpha 1, importin alpha 2, importin subunit alpha-2, importin-alpha-P1, karyopherin subunit alpha-2, pendulin

Official Symbol: KPNA2

Accession Number(s): NP_002257.1

Human GeneID(s): [3838](#)

Non-Human GeneID(s): 16647 (mouse), 85245 (rat)

Immunogen

Peptide with sequence C-QVQDGAPGTFNF., from the C Terminus of the protein sequence according to NP_002257.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Approx 55kDa band observed in lysates of cell line HeLa (calculated MW of 57.9kDa according to NP_002257.1). See non-biotinylated parental product's datasheet for further QC data. Recommended concentration: 0.1-0.3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



Biotinylated EB06233 (0.2 μ g/ml) staining of HeLa lysate (35 μ g protein in RIPA buffer), exactly mirroring its parental non-biotinylated product. Primary incubation was 1 hour. Detected by chemiluminescence, using streptavidin-HRP and using NAP blocker as a substitute for skimmed milk.