



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

## US Office

### Everest Biotech c/o Abcore

405 Maple Street, Suite A106  
Ramona,  
CA 92065  
USA

Inquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[usasales@everestbiotech.com](mailto:usasales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB06237 - Goat Anti-KPNA3 / IPOA4 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** KPNA3, IPOA4, karyopherin alpha 3 (importin alpha 4), RP11-432M24.3, SRP1gamma, SRP4, hSRP1, importin alpha 4, importin alpha 4, importin-alpha-Q2, karyopherin alpha 3, SRP1, importin alpha-3

**Official Symbol:** KPNA3

**Accession Number(s):** NP\_002258.2

**Human GeneID(s):** [3839](#)

**Non-Human GeneID(s):** 16648 (mouse)

### Immunogen

Peptide with sequence C-DPTANLQTKFN, from the C Terminus of the protein sequence according to NP\_002258.2.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:128000.

**Western blot:** Approx. 60kDa band observed in lysates of cell lines CaCo-2, HEK293 and MCF7 (calculated MW of 57.8kDa according to NP\_002258.2). Recommended concentration: 0.01-0.03µg/ml. Primary incubation 1 hour at room temperature.

**IHC:** Paraffin embedded Human Testis. Recommended concentration: 6µg/ml.

**Immunofluorescence:** Strong expression of the protein seen in the cytoplasm and nuclei of A431 cells. Recommended concentration: 10µg/ml.

**Flow Cytometry:** Flow cytometric analysis of A431 cells. Recommended concentration: 10ug/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Mouse, Rat, Dog, Cow

### Specific References

**This antibody has been successfully used in the following paper:**

Sadequl Islam, Yang Sun, Yuan Gao, Tomohisa Nakamura, Arshad Ali Noorani, Tong Li, Philip C Wong, Noriyuki Kimura, Etsuro Matsubara, Kensaku Kasuga, Takeshi Ikeuchi, Taisuke Tomita, Kun Zou, Makoto Michikawa  
Presenilin Is Essential for ApoE Secretion, a Novel Role of Presenilin Involved in Alzheimer's Disease Pathogenesis.  
J Neurosci. 2022 Feb 23;42(8):1574-1586.  
PMID: 34987110

**This antibody (previous batch) has been successfully used in the following paper:**

Krzysztof Sikorski, Adi Mehta, Marit Inngjerdingen, Flourina Thakor, Simon Kling, Tomas Kalina, Tuula A. Nyman, Maria Ekman Stensland, Wei Zhou, Gustavo A. De Souza, Lars

Holden, Jan Stuchly, Markus Templin and Fridtjof Lund-Johansen  
A high-throughput pipeline for validation of antibodies  
Nat Methods. 2018 Nov;15(11):909-912  
PMID: 30377371

**This antibody (previous batch) has been successfully used on Mouse:**


Ana Jovišević, Jerome Mertens, Steven Boeynaems, Elke Bogaert, Noori Chai, Shizuka B. Yamada, Joseph W. Paul III, Shuying Sun, Joseph R. Herdy, Gregor Bieri, Nicholas J. Kramer, Fred H. Gage, Ludo Van Den Bosch, Wim Robberecht, and Aaron D. Gitler.  
Modifiers of C9orf72 dipeptide repeat toxicity connect nucleocytoplasmic transport defects to FTD/ALS.  
Nat Neurosci. 2015 Sep;18(9).  
PMID: 26308983

**This antibody (previous batch) has been successfully used in ICC on Rat:**


Ohara R, Hata K, Yasuhara N, Mehmood R, Yoneda Y, Nakagawa M, Yamashita T.  
Axotomy induces axonogenesis in hippocampal neurons by a mechanism dependent on importin  $\beta$ .  
Biochem Biophys Res Commun. 2011 Feb 25;405(4):697-702.  
PMID: 21291862

**This antibody (previous batch) has been successfully used in WB on Human and Mouse, and in IP on Mouse:**


Huenniger K, Krämer A, Soom M, Chang I, Köhler M, Depping R, Kehlenbach RH, Kaether C.  
Notch1 signaling is mediated by importins alpha 3, 4, and 7.  
Cell Mol Life Sci. 2010 Sep;67(18):3187-96.  
PMID: 20454918




EB06237 (0.03 $\mu$ g/ml) staining of CaCo-2 (A) and HEK293 (B), and (0.01 $\mu$ g/ml) MCF7 (C), cell lysate (35 $\mu$ g protein in RIPA buffer). Detected by chemiluminescence.



EB06237 (6 $\mu$ g/ml) staining of paraffin embedded Human Testis. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.




EB06237 Negative Control showing staining of paraffin embedded Human Testis, with no primary antibody.



EB06237 Immunofluorescence analysis of paraformaldehyde fixed A431 cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing cytoplasmic and nuclear staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).





EB06237 Flow cytometric analysis of paraformaldehyde fixed A431 cells (blue line), permeabilized with 0.5% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (1ug/ml). IgG control: Unimmunized goat IgG (black line) followed by Alexa Fluor 488 secondary antibody.