

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05039 - Goat Anti-LAF4 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: AFF3, AF4/FMR2 family, member 3, LAF4, MLLT2-like, lymphoid nuclear protein related to AF4, MLLT2-related protein, lymphoid nuclear protein 4

Official Symbol: AFF3

Accession Number(s): NP_002276.2; NP_001020279.1

Human GeneID(s): [3899](#)

Immunogen

Peptide with sequence C-SQQGLHWLRNSAHS, from the C Terminus of the protein sequence according to NP_002276.2; NP_001020279.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Western blot: Band observed in Raji at approx 75kDa with 1 hour incubation of primary.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Dog