



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05039 - Goat Anti-LAF4 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** AFF3, AF4/FMR2 family, member 3, LAF4, MLLT2-like, lymphoid nuclear protein related to AF4, MLLT2-related protein, lymphoid nuclear protein 4

**Official Symbol:** AFF3

**Accession Number(s):** NP\_002276.2; NP\_001020279.1

**Human GeneID(s):** [3899](#)

### Immunogen

Peptide with sequence C-SQQGLHWLRNSAHS, from the C Terminus of the protein sequence according to NP\_002276.2; NP\_001020279.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:8000.

**Western blot:** Western blot: Band observed in Raji at approx 75kDa with 1 hour incubation of primary.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Dog