

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06418-T - Goat Anti-LAT1 Antibody - Trial

Size: 20µg specific antibody in 40µl



Target Protein

Principal Names: LAT1, LAT, linker for activation of T cells, pp36, 36 kDa phospho-tyrosine adaptor protein, linker for activation of T-cells family member 1

Official Symbol: LAT

Accession Number(s): NP_055202.1; NP_001014987.1; NP_001014989.1;

NP_001014988.1

Human GenelD(s): 27040

Important Comments: This antibody is expected to recognize isoform a (NP_055202.1), isoform b (NP_001014987.1 and NP_001014989.1) and isoform c (NP_001014988.1).

Immunogen

Peptide with sequence C-GAPDYENLQELN, from the C Terminus of the protein sequence according to NP_055202.1; NP_001014987.1; NP_001014989.1; NP_001014988.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:2000.

Western blot: Approx 37kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 27.9kDa according to NP_055202.1). This molecular weight is routinely observed by other sources. Recommended concentration: 0.01-0.03μg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat



EB06418 (0.01µg/ml) staining of Peripheral Blood Mononucleocytes lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.