



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08116 - Goat Anti-Leucine-rich repeat kinase 1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: LRRK1, leucine-rich repeat kinase 1, FLJ23119, FLJ27465, KIAA1790, RIPK6, Roco1

Official Symbol: LRRK1

Accession Number(s): NP_078928.3

Human GeneID(s): [79705](#)

Non-Human GeneID(s): 233328 (mouse)

Immunogen

Peptide with sequence C-ELPTEPTDDNP, from the internal region of the protein sequence according to NP_078928.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Western Blot: Preliminary experiments gave an approx 38kDa band in Human Brain (Cerebral Cortex, Frontal Cortex) and Prostate lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 225kDa according to NP_078928.3. The 38kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat,