

### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB08099 - Goat Anti-LITAF Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

**Principal Names:** LITAF, lipopolysaccharide-induced TNF factor, CMT1C, FLJ38636, MGC116698, MGC116700, MGC116701, MGC125274, MGC125275, MGC125276, PIG7, SIMPLE, TP53I7, LPS-induced TNF-alpha factor, small integral membrane protein of lysosome/late endosome, tumor protein

Official Symbol: LITAF

Accession Number(s): NP\_004853.2; NP\_001129945.1

Human GeneID(s): 9516

Non-Human GenelD(s): 56722 (mouse), 65161 (rat)

Important Comments: This antibody is excpected to crossreact with both reported

isoforms (NP\_004853.2; NP\_001129945.1).

# **Immunogen**

Peptide with sequence C-PDGKGMNPPSYYTQ, from the internal region of the protein sequence according to NP\_004853.2; NP\_001129945.1.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:8000.


**Western blot:** Preliminary experiments gave an approx 25kDa band in Human Lymph Node lysates and in lysates of cell lines HeLa and HepG2 after 0.1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 17.1kDa according to NP\_004853.2. The 25kDa band was successfully blocked by incubation with the immunizing peptide.

IHC: Paraffin embedded Human Thymus. Recommended concentration: 6.5µg/ml.

# **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat



EB08099 (6.5μg/ml) staining of paraffin embedded Human Thymus. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.