

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06808 - Goat Anti-Mammaglobin B / SCGB2A1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: SCGB2A1, secretoglobin, family 2A, member 1, HGNC:7051, LPHC,

MGB2, MGC71973, UGB3, lipophilin C, mammaglobin 2, mammaglobin B

Official Symbol: SCGB2A1

Accession Number(s): NP_002398.1

Human GeneID(s): 4246

Immunogen

Peptide with sequence CFLNQSHRTLKN, from the internal region of the protein sequence according to NP_002398.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:2000.

Species Reactivity

Tested:

Expected from sequence similarity: Human