



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB08302 - Goat Anti-MARCH-VII / Axotrophin Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MARCH7, membrane-associated ring finger (C3HC4) 7, AXO, AXOT, DKFZP586F1122, MARCH-VII, RNF177, axotrophin, membrane-associated RING-CH protein VII

Official Symbol: MARCH7

Accession Number(s): NP_073737.1

Human GeneID(s): [64844](#)

Non-Human GeneID(s): 57438 (mouse), 311059 (rat)

Immunogen

Peptide with sequence C-NFNQESEGRNT, from the internal region of the protein sequence according to NP_073737.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.


Western blot: Preliminary experiments gave an approx 125kDa band in Human Spleen and Tonsil lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 78.1kDa according to NP_073737.1. The 125kDa band was successfully blocked by incubation with the immunizing peptide.

IHC: Paraffin embedded Human Uterus. Recommended concentration: 3.75µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat



EB08302 (3.75µg/ml) staining of paraffin embedded Human Uterus. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.