



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05252 - Goat Anti-MPP2 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** MPP2, membrane protein, palmitoylated 2 (MAGUK p55 subfamily member 2), DLG2, DKFZp761D0712, discs large, homolog 2, MAGUK p55 subfamily member 2

**Official Symbol:** MPP2

**Accession Number(s):** NP\_005365

**Human GeneID(s):** [4355](#)

### Immunogen

Peptide with sequence PVAATNSETAMQQ-C, from the N Terminus of the protein sequence according to NP\_005365.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Western blot:** Western blot: Two consistent bands observed in HeLa at approx 70kDa and approx 50kDa with 12 hours incubation of primary.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human