



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05252 - Goat Anti-MPP2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MPP2, membrane protein, palmitoylated 2 (MAGUK p55 subfamily member 2), DLG2, DKFZp761D0712, discs large, homolog 2, MAGUK p55 subfamily member 2

Official Symbol: MPP2

Accession Number(s): NP_005365

Human GeneID(s): [4355](#)

Immunogen

Peptide with sequence PVAATNSETAMQQ-C, from the N Terminus of the protein sequence according to NP_005365.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Two consistent bands observed in HeLa at approx 70kDa and approx 50kDa with 12 hours incubation of primary.

Species Reactivity

Tested:

Expected from sequence similarity: Human