

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB09331 - Goat Anti-MRGPRX3 / MRGX3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MRGPRX3, MAS-related GPR, member X3, GPCR, MRGX3, SNSR1,

G protein-coupled receptor MRGX3, G protein-coupled receptor SNSR1, G

protein-coupled receptor SNSR2, mas-related GPCR member X3

Official Symbol: MRGPRX3

Accession Number(s): NP_473372.2

Human GeneID(s): 117195

Important Comments: This antibody may cross-react with human MRGPRX1 and X4

(GeneID: 259249, 117196; MW: 36.25kDa and 36.43kDa).

Immunogen

Peptide with sequence C-QLPQETLELSGSRL, from the Internal region of the protein sequence according to NP_473372.2.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Species Reactivity

Tested:

Expected from sequence similarity: Human