



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05008 - Goat Anti-MRVI1 / IRAG Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MRVI1, murine retrovirus integration site 1 homolog, IRAG, inositol 1,4,5-triphosphate-associated cGMP kinase substrate

Official Symbol: MRVI1

Accession Number(s): NP_569056; NP_006060

Human GeneID(s): [10335](#)

Important Comments: This antibody is expected to recognise both reported isoforms (NP_569056; NP_006060).

Immunogen

Peptide with sequence C-WSSGLQHEQPTEQ, from the C Terminus of the protein sequence according to NP_569056; NP_006060.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Species Reactivity

Tested:

Expected from sequence similarity: Human