



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05164 - Goat Anti-MTA2 / MTA1L1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MTA1L1, metastasis-associated 1-like 1, PID, MTA2, MTA1-L1, metastasis-associated protein 2, metastasis associated gene 1-like 1, DKFZp686F2281, metastasis associated gene family, member 2, metastasis-associated 1-like 1

Official Symbol: MTA1L1

Accession Number(s): NP_004730.2

Human GeneID(s): [9219](#)

Immunogen

Peptide with sequence C-LKVKPTLIAVRP, from the C Terminus of the protein sequence according to NP_004730.2.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Western blot: Approx 75-85 kDa band seen in Jurkat and Hela Lysate. Recommended for use at 2-4 µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse