



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB05238 - Goat Anti-MTMR1 (isoform 1) Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** MTMR1, myotubularin related protein 1, myotubularin-related protein 1

**Official Symbol:** MTMR1

**Accession Number(s):** NP\_003819.1

**Human GeneID(s):** [8776](#)

**Important Comments:** Please note this will only recognise Human isoform 1 (NP\_003819.1), not isoform 2 (NP\_789746).

### Immunogen

Peptide with sequence C-SSPSHSATSVHTSV, from the C Terminus of the protein sequence according to NP\_003819.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human