



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05193 - Goat Anti-MTMR3 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: MTMR3, myotubularin related protein 3, KIAA0371, FYVE-DSP1, FYVE (Fab1 YGLO23 Vsp27 EEA1 domain) dual-specificity protein phosphatase

Official Symbol: MTMR3

Accession Number(s): NP_066576

Human GeneID(s): [8897](#)

Immunogen

Peptide with sequence C-DLELDKPIAATSN, from the C Terminus of the protein sequence according to NP_066576.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: No signal obtained yet but low background observed in Mouse Brain and HeLa lysates at upto 3µg/ml.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse