

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB11103 - Goat Anti-Muc16 (Rat) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: 1110008I14Rik, CA125, CA-125, CA125 ovarian cancer antigen, Gm1480, MUC-16, mucin 16, cell surface associated, mucin-16, ovarian cancer-related tumor marker CA125, ovarian carcinoma antigen CA125, RGD1308422, MUC16

Official Symbol: Muc16

Accession Number(s): XP_235886.5

Non-Human GenelD(s): 73732 (mouse), 315451 (rat)

Immunogen

Peptide with sequence C-QRHRLAYYLSHLDLRK, from the internal region of the protein sequence according to XP_235886.5.

Please note the <u>peptide</u> is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:1000.

Western blot: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 794kDa according to XP_235886.5. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring product. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Rat, Mouse