



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

US Office

Everest Biotech c/o Abcore

405 Maple Street, Suite A106
Ramona,
CA 92065
USA

Inquiries:

info@everestbiotech.com

Sales:

usasales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06140 - Goat Anti-MUNC18 / STXBP1 (isoform a) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: STXBP1, MUNC18-1, UNC18, hUNC18, syntaxin binding protein 1, Hunc18, rbSec1, syntaxin-binding protein 1, RP11-56D16.3, EIEE4, FLJ37475, RBSEC1

Official Symbol: STXBP1

Accession Number(s): NP_003156.1

Human GeneID(s): [6812](#)

Non-Human GeneID(s): 20910 (mouse)

Important Comments: This antibody is expected to recognize isoform a (NP_003156.1).

Immunogen

Peptide with sequence C-SRVSFEDQAPTME, from the C Terminus of the protein sequence according to NP_003156.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Approx 70-75kDa band observed in Human Cerebellum lysates (calculated MW of 68.7kDa according to Human NP_003156.1). Recommended concentration: 0.1-0.3µg/ml. Primary incubation 1 hour at room temperature.

IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse

Specific Reference

This antibody has been successfully used in IHC on Mouse:

Meijer M, Cijssouw T, Toonen RF, Verhage M.


Synaptic Effects of Munc18-1 Alternative Splicing in Excitatory Hippocampal Neurons.

PLoS One. 2015 Sep 25;10(9):e0138950.


PMID: 26407320



EB06140 (0.3 μ g/ml) staining of Human Cerebellum lysate (35 μ g protein in RIPA buffer). Detected by chemiluminescence.



EB06140 (5 μ g/ml) staining of paraffin embedded Human Cortex. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.



EB06140 Negative Control showing staining of paraffin embedded Human Cortex, with no primary antibody.