

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB10793 - Goat Anti-NDRG2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: NDRG2, NDRG family member 2, DKFZp781G1938, FLJ25522, KIAA1248, SYLD, N-myc downstream regulator 2, NDR1-related protein NDR2, OTTHUMP00000164405, OTTHUMP00000164407, OTTHUMP00000164408,

cytoplasmic protein Ndr1, protein NDRG2, syld709

Official Symbol: NDRG2

Accession Number(s): NP_963293.1

Human GeneID(s): 57447

Non-Human GenelD(s): 29811 (mouse), 171114 (rat)

Important Comments: This antibody is expected to recognize reported isoform a (NP_963293.1) only. Reported variants represent identical protein: NP_963831.1,

NP_963293.1, NP_963834.1, NP_963833.1.

Immunogen

Peptide with sequence KEAELAARILLDQGQ-C, from the N Terminus (near) of the protein sequence according to NP_963293.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Approx 45kDa band observed in Human Brain (Temporal Cortex), Mouse and Rat Brain lysates (calculated MW of 40.8kDa according to NP_963293.1).


Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.

IHC: Paraffin embedded Human Adrenal Gland. Recommended concentration: 3.75µg/ml.


Species Reactivity

Tested: Human, Mouse, Rat


Expected from sequence similarity: Human, Mouse, Rat, Dog



EB10793 (0.1µg/ml) staining of Human Temporal Cortex lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.



EB10793 (0.3µg/ml) staining of Mouse (A) and Rat (B) Brain lysates (35µg protein in RIPA buffer). Detected by chemiluminescence.



EB10793 (3.75µg/ml) staining of paraffin embedded Human Adrenal Gland. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.