



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05085 - Goat Anti-NESCA Adapter protein Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: RUSC1, NESCA, RUN and SH3 domain containing 1, nesca protein, new molecule containing SH3 at the carboxy-terminus

Official Symbol: RUSC1

Accession Number(s): NP_055143

Human GeneID(s): [23623](#)

Immunogen

Peptide with sequence C-EGLVPVGYTSLVL, from the C Terminus of the protein sequence according to NP_055143.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Western blot: Western blot: Approx 50kDa band seen in HeLa and U937 cell lysates. Recommended for use at 1-3 µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human