



## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

Fax: +44 (0)1869 238327

## US Office

### Everest Biotech c/o Abcore

405 Maple Street, Suite A106  
Ramona,  
CA 92065  
USA

Inquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[usasales@everestbiotech.com](mailto:usasales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: 888-320-4628 (toll-free)

Fax: 888-841-9041

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB08482 - Goat Anti-Neurotrophin 5 / NT4 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** NT4, neurotrophin 5 (neurotrophin 4/5), NTF5, NT-4/5, NT5, NTF4, neurotrophic factor 4, neurotrophic factor 5, neurotrophin 4, neurotrophin 5

**Official Symbol:** NTF5

**Accession Number(s):** NP\_006170.1

**Human GeneID(s):** [4909](#)

### Immunogen

Peptide with sequence FETRCKADNAEEG, from the internal region of the protein sequence according to NP\_006170.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested


**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** Approx 55kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 22.4kDa according to NP\_006170.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Ghinelli et al, Invest Ophthalmol Vis Sci. 2003 Aug;44(8):3352-7; PMID: 12882781). Recommended concentration: 0.1-0.3µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human



EB08482 (0.1µg/ml) staining of Human Cerebellum lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.