

## UK Office

### Everest Biotech Ltd

Cherwell Innovation Centre  
77 Heyford Park  
Upper Heyford  
Oxfordshire  
OX25 5HD  
UK

Enquiries:

[info@everestbiotech.com](mailto:info@everestbiotech.com)

Sales:

[sales@everestbiotech.com](mailto:sales@everestbiotech.com)

Tech support:

[support@everestbiotech.com](mailto:support@everestbiotech.com)

Tel: +44 (0)1869 238326

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB08176 - Goat Anti-NKD2 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** NKD2, naked cuticle homolog 2 (Drosophila), Dvl-binding protein  
NKD2, naked cuticle homolog 2, naked cuticle-2

**Official Symbol:** NKD2

**Accession Number(s):** NP\_149111.1

**Human GeneID(s):** [85409](#)

**Non-Human GeneID(s):** 72293 (mouse), 308068 (rat)

### Immunogen

Peptide with sequence C-HKRYRQKGREGHS, from the internal region of the protein  
sequence according to NP\_149111.1.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity  
chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum  
albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

### Species Reactivity

**Tested:**

**Expected from sequence similarity:** Human, Mouse, Rat