



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06758 - Goat Anti-Nogo 66 receptor Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Nogo-66 receptor, RTN4R, reticulon 4 receptor, HGNC:1860, NGR, NOGOR, UNQ330/PRO526, nogo receptor, reticulon 4 receptor precursor, OTTHUMP00000197866

Official Symbol: RTN4R

Accession Number(s): NP_075380.1

Human GeneID(s): [65078](#)

Immunogen

Peptide with sequence C-LQALPDDTFRDLGN, from the internal region of the protein sequence according to NP_075380.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Experiments gave bands at approx 90-100kDa and 37+40+45kDa in Human brain lysate after 0.3µg/ml antibody staining. These bands were successfully blocked by the immunizing peptide and they correspond to earlier findings in the literature with different antibodies (Wang et al, J Neurosci. 2002 Jul 1;22(13):5505-15). This protein has a calculated MW of 50.7kDa according to NP_075380.1. Recommended concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Pig, Dog



EB06758 (0.3µg/ml) staining of Human Brain lysate (35µg protein in RIPA buffer) with (B) and without (A) blocking with the immunising peptide. Primary incubation was 1 hour. Detected by chemiluminescence.