

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06758 - Goat Anti-Nogo 66 receptor Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Nogo-66 receptor, RTN4R, reticulon 4 receptor, HGNC:1860, NGR,

NOGOR, UNQ330/PRO526, nogo receptor, reticulon 4 receptor precursor,

OTTHUMP00000197866
Official Symbol: RTN4R

Accession Number(s): NP_075380.1

Human GeneID(s): 65078

Immunogen

Peptide with sequence C-LQALPDDTFRDLGN, from the internal region of the protein sequence according to NP_075380.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.


Western blot: Experiments gave bands at approx 90-100kDa and 37+40+45kDa in Human brain lysate after 0.3μg/ml antibody staining. These bands were successfully blocked by the immunizing peptide and they correspond to earlier findings in the literature with different antibodies (Wang et al, J Neurosci. 2002 Jul 1;22(13):5505-15). This protein has a calculated MW of 50.7kDa according to NP_075380.1. Recommended

concentration: 0.3-1µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Pig, Dog



EB06758 (0.3 μ g/ml) staining of Human Brain lysate (35 μ g protein in RIPA buffer) with (B) and without (A) blocking with the immunising peptide. Primary incubation was 1 hour. Detected by chemiluminescence.