

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10096 - Goat Anti-Npc1(mouse) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: A430089E03Rik, C85354, D18Ertd139e, D18Ertd723e, lcsd, Niemann Pick type C1, niemann-Pick C1 protein, nmf164, sphingomyelinosis, spm, Npc1

Official Symbol: Npc1

Accession Number(s): NP_032746.2

Non-Human GeneID(s): 18145 (mouse)

Immunogen

Peptide with sequence C-TYHTILKTSADYTD, from the internal region of the protein sequence according to NP_032746.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Preliminary experiments gave an approx 36kDa band in lysates of cell lines NIH3T3 and NSO after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 142kDa according to NP_032746.2. The 36kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?

Species Reactivity

Tested:

Expected from sequence similarity: Mouse