



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10837 - Goat Anti-Npc111 (mouse) (aa56-69) Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: Npc111, NPC1-like 1, RP23-198N14.2, 9130221N23Rik, Gm243, Niemann-Pick C1-like protein 1, Niemann-Pick disease, type C1, Niemann-Pick disease, type C1-like 1, OTTMUSP00000000725

Official Symbol: Npc111

Accession Number(s): NP_997125.2

Non-Human GeneID(s): 237636 (mouse), 432367 (rat)

Immunogen

Peptide with sequence CLSNTPARHVTGDH, from the N Terminus of the protein sequence according to NP_997125.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Preliminary experiments in Mouse and Rat Liver lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Species Reactivity

Tested:

Expected from sequence similarity: Mouse, Rat