



International Office

Everest Biotech Ltd

Vector Laboratories, Inc.
6737 Mowry Ave
Newark, CA 94560
United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB10400 - Goat Anti-NRAP Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: nebulin-related anchoring protein, NRAP

Official Symbol: NRAP

Accession Number(s): NP_006166.3; NP_932326.2

Human GeneID(s): [4892](#)

Non-Human GeneID(s): 18175 (mouse), 307982 (rat)

Important Comments: This antibody is expected to recognize all reported isoforms (NP_006166.3; NP_932326.2).

Immunogen

Peptide with sequence EDTKSQCHVSLD, from the internal region of the protein sequence according to NP_006166.3; NP_932326.2.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:4000.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow, Pig