



## International Office

### Everest Biotech Ltd

Vector Laboratories, Inc.  
6737 Mowry Ave  
Newark, CA 94560  
United States

Customer Service:

[customerservice@vectorlabs.com](mailto:customerservice@vectorlabs.com)

Technical Service:

[technical@vectorlabs.com](mailto:technical@vectorlabs.com)

Tel: +1 (800) 227-6666

[www.everestbiotech.com](http://www.everestbiotech.com)

**Research Use Only. Not for  
diagnostic or therapeutic use.**

## EB08418 - Goat Anti-OCT6 Antibody

Size: 100µg specific antibody in 200µl



### Target Protein

**Principal Names:** POU3F1, POU class 3 homeobox 1, OCT6, OTF6, SCIP, POU domain, class 3, transcription factor 1, octamer-binding transcription factor 6

**Official Symbol:** POU3F1

**Accession Number(s):** NP\_002690.3

**Human GeneID(s):** [5453](#)

### Immunogen

Peptide with sequence C-HHALHEDGHEAQLE, from the internal region of the protein sequence according to NP\_002690.3.

Please note the [peptide](#) is available for sale.

### Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

### Applications Tested

**Peptide ELISA:** antibody detection limit dilution 1:32000.

**Western blot:** Approx 50kDa band observed in Human Ovary, Placenta and Umbilical Cord lysates (calculated MW of 45.5kDa according to NP\_002690.3). Recommended concentration: 1-3µg/ml.

### Species Reactivity

**Tested:** Human

**Expected from sequence similarity:** Human, Mouse, Rat, Dog, Cow



EB08418 (2µg/ml) staining of Human Umbilical Cord lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.