

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB08418 - Goat Anti-OCT6 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

**Principal Names:** POU3F1, POU class 3 homeobox 1, OCT6, OTF6, SCIP, POU domain, class 3, transcription factor 1, octamer-binding transcription factor 6

Official Symbol: POU3F1

Accession Number(s): NP\_002690.3

Human GeneID(s): 5453

### **Immunogen**

Peptide with sequence C-HHALHEDGHEAQLE, from the internal region of the protein sequence according to NP\_002690.3.

Please note the peptide is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Aliquot and store at -20°C. Minimize freezing and thawing.

#### **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:32000.

**Western blot:** Approx 50kDa band observed in Human Ovary, Placenta and Umbilical Cord lysates (calculated MW of 45.5kDa according to NP\_002690.3). Recommended

concentration: 1-3µg/ml.

### **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow



EB08418 (2μg/ml) staining of Human Umbilical Cord lysate (35μg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.