



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06207 - Goat Anti-OCT6 / SCIP / POU3F1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: POU3F1, OCT6, OTF6, OCT-6, POU domain, class 3, transcription factor 1, SCIP, POU class 3 homeobox 1, octamer-binding transcription factor 6

Official Symbol: POU3F1

Accession Number(s): NP_002690.3

Human GeneID(s): [5453](#)

Non-Human GeneID(s): 18991 (mouse), 192110 (rat)

Immunogen

Peptide with sequence ATTAQYLPRGPGG-C, from the N Terminus of the protein sequence according to NP_002690.3.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Species Reactivity

Tested:

Expected from sequence similarity: Human, Mouse, Rat, Cow