

International Office

Everest Biotech Ltd

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

Customer Service:

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB09023 - Goat Anti-OXT Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: OXT, oxytocin, prepro- (neurophysin I), MGC126890, MGC126892,

OT, OT-NPI, oxytocin-neurophysin I, oxytocin-neurophysin I, preproprotein

Official Symbol: OXT

Accession Number(s): NP_000906.1

Human GeneID(s): 5020

Immunogen

Peptide with sequence HADPACDAEATFSQR, from the C Terminus of the protein sequence according to NP_000906.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:4000.

IHC: Frozen section of Human Hypothalamus shows staining of magnocellular oxytocin neurons in the supraoptic nucleus. Recommended concentration, 0.05-0.15µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB09023 (0.05µg/ml) staining of PFA-perfused cryosection of Human Hypothalamus. Antigen retrieval with citrate buffer pH 6 at 80C for 30min, HRP-staining with Ni-DAB after Biotin-SP-antigoat amplification. Data obtained by Prof. Erik Hrabovszky, Inst, Exp, Med., Budapest, Hungary.