

### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB05630 - Goat Anti-p70S6K / RPS6KB1 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

**Principal Names:** RPS6KB1, p70-S6K, ribosomal protein S6 kinase, 70kD, polypeptide 1, PS6K, S6K, S6K1, STK14A, p70(S6K)-alpha, p70-alpha, p70 S6 kinase, alpha 1, p70 S6 kinase, alpha 2, ribosomal protein S6 kinase, 70kD, polypeptide 1, serine/threonine kinase

14 alpha

Official Symbol: RPS6KB1

Accession Number(s): NP\_003152.1

Human GeneID(s): 6198

Non-Human GenelD(s): 72508 (mouse), 83840 (rat)

## Immunogen

Peptide with sequence C-MISKRPEHLRMNL, from the C Terminus of the protein sequence according to NP\_003152.1.

Please note the peptide is available for sale.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:64000.

**Western blot:** Approx 60kDa band observed in lysates of Human Placenta and of cell line NIH3T3 (calculated MW of 59.1kDa according to NP\_003152.1). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour. Preliminary testing was unsuccessful on Rat Kidney for this particular batch.


#### **Species Reactivity**

Tested: Human, Mouse

Expected from sequence similarity: Human, Mouse, Rat, Dog, Cow



EB05630 (0.3 $\mu$ g/ml) staining of Human Placenta lysate (35 $\mu$ g protein in RIPA buffer). Detected by chemiluminescence.



 $EB05630\ (1\mu\text{g/ml})\ staining\ of\ NIH3T3\ cell\ lysate\ (35\mu\text{g}\ protein\ in\ RIPA\ buffer).\ Detected\ by\ chemiluminescence.$