

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB06284 - Goat Anti-P8 / NUPR1 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

**Principal Names:** NUPR1, Nuclear Protein 1, P8, COM1, p8 protein (candidate of metastasis 1), p8 protein, nuclear protein, transcriptional regulator, 1, candidate of

metastasis 1

Official Symbol: NUPR1

Accession Number(s): NP\_001035948.1; NP\_036517

Human GeneID(s): 26471

### **Immunogen**

Peptide with sequence ATFPPATSAPQQP-C, from the N Terminus of the protein sequence according to NP\_001035948.1; NP\_036517.

Please note the <u>peptide</u> is available for sale.

# **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

Peptide ELISA: antibody detection limit dilution 1:32000.

### **Species Reactivity**

Tested:

Expected from sequence similarity: Human