

#### **International Office**

#### **Everest Biotech Ltd**

Vector Laboratories, Inc. 6737 Mowry Ave Newark, CA 94560 United States

**Customer Service:** 

customerservice@vectorlabs.com

Technical Service:

technical@vectorlabs.com

Tel: +1 (800) 227-6666

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

# EB05531 - Goat Anti-PAIP2 Antibody

Size: 100µg specific antibody in 200µl



## **Target Protein**

Principal Names: PAIP2, PABP-interacting protein 2, polyadenylate binding

protein-interacting protein

Official Symbol: PAIP2

Accession Number(s): NP\_057564

Human GeneID(s): 51247

Non-Human GenelD(s): 67869 (mouse)

#### **Immunogen**

Peptide with sequence KDPSRSSTSPSIIN-C, from the N Terminus of the protein sequence according to NP\_057564.

Please note the peptide is available for sale.

## **Purification and Storage**

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

## **Applications Tested**

**Western blot:** Western blot: Approx 15kDa band seen in 293 lysate. Recommended for use at 1-3µg/ml.

## **Species Reactivity**

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat