



UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09719 - Goat Anti-PCDH11Y / PCDH11X Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PCDH11Y, protocadherin 11 Y-linked, PCDH-PC, PCDH11X, PCDH22, PCDHY, OTTHUMP00000033077, protocadherin 11 X-linked, protocadherin 22, protocadherin Y, RP13-138P15.1, PCDH-X, PCDH11, PCDHX, protocadherin 11X, protocadherin X, protocadherin-S

Official Symbol: PCDH11Y

Accession Number(s): NP_116753.1; NP_116754.1; NP_116755.1; NP_055337.1; NP_055337.1; NP_116750.1; NP_116751.1

Human GeneID(s): [83259](#) , [27328](#)

Important Comments: This antibody is expected to recognize all reported Y-linked isoforms (NP_116753.1; NP_116754.1; NP_116755.1) and all reported X-linked isoforms (NP_055337.1; NP_055337.1; NP_116750.1; NP_116751.1).

Immunogen

Peptide with sequence C-KYNWVTTPTTFK, from the internal region of the protein sequence according to NP_116753.1; NP_116754.1; NP_116755.1; NP_055337.1; NP_055337.1; NP_116750.1; NP_116751.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.


Western blot: Preliminary experiments gave an approx 75kDa band in Human Brain (Amygdala, Hippocampus, Substantia nigra) lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 115kDa according to NP_116753.1 and of 113kDa according to NP_055337.1, the smallest reported isoforms. The 75kDa band was successfully blocked by incubation with the immunizing peptide.

IHC: Paraffin embedded Human Brain (Cortex). Recommended concentration: 6µg/ml.


Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB09719 (6 μ g/ml) staining of paraffin embedded Human Cortex. Heat induced antigen retrieval with citrate buffer pH 6, HRP-staining.



EB09719 Negative Control showing staining of paraffin embedded Human Cortex, with no primary antibody.